# Package: tgram (via r-universe)

October 18, 2024

Date/Publication 2017-11-24 23:31:50 UTC

# Contents

Index

juniperus	2	2
standz.all	2	2
tgram	4	ŀ
traq.profile	6	5
	7	1

juniperus

#### Description

An example of traqueid measurements to standarize with function tgram.

#### Usage

data(juniperus)

#### Format

A data frame with 77 observations on the following 4 variables.

traqueidogram Numeric vector indicating the traqueidogram to which each measurement belongs lumen.wall A factor indicating if the measurement is lumen (1)or wall (w) order Position of the measurement in the ordered sequence within each traqueidogram

width.um Width (micrometres) of each measurement

## Examples

standz.all Vaganov Normalized Tracheidogram

#### Description

The function produces a normalized tracheidogram, i.e. a curve showing variations in cell parameters as a function of the cell position within an annual ring, following the procedure of Vaganov (1990).

#### Usage

```
standz.all(traq, series, wl = NULL, w.char = NULL, order = NULL, G = 30)
standz(tgl1, G=30)
## S3 method for class 'standz.all'
plot(x, which=NULL,...)
```

# standz.all

#### Arguments

traq	A vector with the ordered sequences of measurements for each tracheidogram.
series	A vector of indicator values (i.e. a factor) with each level indicating each unique tracheidogram.
wl	A vector indicating if the measurement is wall or lumen.
w.char	Character used in w1 to indicate "wall".
order	Vector indicating the ordering of each measurement in each lumen or wall series within a tracheidogram.
G	Number of cells to get the original measurements normalized to.
tgl1	Vector with the ordered sequences of measurements of a single tracheidogram.
х	An object of class standz.all resulting from applying the standz.all function.
which	One of NULL, "w", or any other character. This produces the plot function to draw all the tracheidograms together, only the "wall" traqueidograms or only the "lumen" ones, respectively.
	Additional graphical parameters passed to link{plot}.

# Value

standz returns a vector of length G with the normalized values. standz.all returns an object of class standz.all. Basically a list with the following elements:

data.stdz	A matrix with <i>G</i> columns and as many rows as unique wall and lumen tra- cheidograms were in the original data, each with the normalized values of each tracheidogram.
which.w	Vector indicating which rows in data.stdz are "wall" tracheidograms.
which.l	Vector indicating which rows in data.stdz are "lumen" tracheidograms.

# Author(s)

Marcelino de la Cruz Rot and Lucia DeSoto

#### References

Vaganov, E.A. 1990. The tracheidogram method in tree-ring analysis and its application. In: Cook E.R., Kairiukstis L.A., eds. *Methods of dendrochronology: applications in the environmental sciences*. Kluwer Academic Publishers. Dordrecht, the Netherlands. pp. 63-76.

# Examples

```
plot(cosa, type="1")
plot(cosa, type="1", which="w")
plot(cosa, type="1", which="1", add=TRUE)
standz(with(juniperus,width.um[traqueidogram==1 & lumen.wall =="1"]), G=20)
lines(1:20,
    standz(with(juniperus,width.um[traqueidogram==1 & lumen.wall =="1"]), G=20),
    lwd=3)
```

tgram

Compute Tracheidograms

# Description

Function to compute tracheidograms from microscopic measurements.

# Usage

```
tgram(traq, val50 = 50, xlim = NULL, ylim = NULL, mw = 1, plotit = TRUE)
```

```
## S3 method for class 'tgram'
plot(x, xlim = NULL, ylim = NULL, colores = c("red", "green"),
    leyenda = c("lumen", "double wall"), lwd = 2,
    add = FALSE,traq.0 = TRUE, bg.legend = NULL, ...)
```

#### Arguments

traq	microscopic light profile. A vector with the light measurements or a data.frame or matrix with two columns. In this case, the first column contains pixel order (i.e. position along the x axis) and the second one contains the light (i.e. "y") lectures.
val50	"y-value" at wich measurements should be made.
mw	integer. Width of the rolling window to smooth the original data.
plotit	logical. should the original data and the computed tracheidogram be drawn?
xlim	extent of the x-axis. A vector of length 2.
ylim	extent of the y-axis. A vector of length 2.
х	an object of class tgram, resulting from tgram function.
colores	a vector of length 2, with the colors to draw the lumen and wall measurements, respectively.
leyenda	a vector of length 2 with the legend to appear in the plot. By default leyenda =c("lumen","double wall").
lwd	width of the lines in the legend.
add	logical. If TRUE, add to a current plot.

4

traq.0	logical. If TRUE, draw the original measurements.
bg.legend	background color for the legend.
	additional graphical parameters.

#### Details

The purpouse of this function is obtaining cell anatomical data from microscopic light measurements (see DeSoto et al. for details of data adquisition). The microscopic lectures are first smoothed with a rolling window (using function rollmean of package zoo, using the selected width mw. Then, the smoothed curve is "cut" at the treshold value val50 and the distances among the intersection points are computed. This provides an ordered sequence of lumen diameters (LD) and double wall thikness (DWT) measurements. From this sequence some other anatomical measurements are computed. Radial cell wall thickness is computed as CWT[t] = 1/2 \* (DWT[t]/2 + DWT[t+1]/2). Tracheid diameter is computed as TD[t] = DWT[t]/2 + LD[t] + DWT[t+1]/2.

### Value

tgram returns an object of class tgram, bassically a list with

traq	original data.
traq0	if traq was a 2 column matrix, then traq0 returns the same object. If traq was a vector, traq0 returns a two column matrix (first column with pixel position and second with ligh lectures).
cut.points	two column matrix with the coordinates of the intersection of $y = val50$ and the smoothed curve.
what	vector indicating if the computed distances are of lumen $(1)$ or double wall $(2)$ .
distances	ordered sequence of the computed distances (both of lumen and double wall).
LD	ordered sequence of lumen diameters.
DWT	ordered sequence of double wall thickness.
mw	width of the rolling window employed to smooth the data.
CWT	ordered sequence of radial cell wall thikness'.
TD	ordered sequence of tracheid diameters.
LD_CWT_ratio	ordered sequence of LD/CWT ratio.

#### Author(s)

Marcelino de la Cruz Rot and Lucia DeSoto

#### References

DeSoto, L., De la Cruz, M. & Fonti, P. 2011. Intra-annual pattern of tracheid size in the Mediterranean Juniperus thurifera as indicator for seasonal water stress. *Canadian Journal of Forest Research* 41: 1280-1294.

# Examples

data(traq.profile)

traq.profile

Light Throughout a Microscopic Section of Juniperus Wood

# Description

A measurement of light intensity troughout a sequence of pixels in a microscopic section of *Juniperus thurifera* wood.

# Usage

data(traq.profile)

#### Format

A data frame with 883 observations on the following 2 variables.

X1 Pixel position in the sequence

Y1 Light intensity

# Examples

```
data(traq.profile)
tgram(traq.profile, mw=10)
```

6

# Index

\* datasets juniperus, 2 traq.profile,6 \* manip standz.all,2 tgram,4 \* smooth standz.all, 2tgram,<mark>4</mark> juniperus, 2 plot.standz.all(standz.all), 2 plot.tgram(tgram), 4 rollmean, 5 standz(standz.all), 2 standz.all, 2 tgram,<mark>4</mark> traq.profile, 6